

Western Blotting

After determining cell lysate concentration, lysates (total of 50 ug protein) were mixed with sample buffer (two volumes) and heated on the heat block at 90 C for 10 min. Centrifuged, put on ice and loaded on gel.

SDS-PAGE Running Buffer (Towbin)- 2 L

25 mM Tris, 192 mM glycine, 0.1% SDS

1X Running Buffer

Reagents needed:

28.8 g glycine

6.04 g Tris base

2 g SDS

1.8 L ddH₂O

10X Running Buffer

Reagents needed:

288 g glycine

60.4 g Tris base

20 g SDS

1.8 L ddH₂O

**** CAUTION ** SDS powder is hazardous. Prepare solution in a ventilated fume hood.**

Directions:

- 1) Dissolve Tris base and glycine together in 1.8 L of ddH₂O.
- 2) Add SDS and mix.
- 3) Add ddH₂O to a final volume of 2 L.

Transfer buffer.

25 mM Tris, 192 mM glycine, 10% methanol

1X Transfer Buffer

Reagents needed:

28.8 g glycine

6.04 g Tris base

200 ml methanol

1.6 L ddH₂O

10X Transfer Buffer

Reagents needed:

288 g glycine

60.4 g Tris base

- methanol

1.8 L ddH₂O

**** NOTE: for the proper transfer of large proteins, up to 0.5% SDS may need to be added to 1X Transfer Buffer. ****

Directions for 1X Transfer Buffer:

- 1) Dissolve Tris base and glycine together in 1.6 L of ddH₂O.
- 2) Add methanol and mix.
- 3) Add ddH₂O to a final volume of 2 L.

Directions for 10X Transfer Buffer:

- 1) Dissolve Tris base and glycine together in 1.8 L of ddH₂O.
- 2) Add ddH₂O to a final volume of 2 L.

**** To make 1X Transfer Buffer from 10X: Mix 100 ml of 10X Transfer Buffer, 100 ml of methanol and 800 ml of ddH₂O per liter ****

Membrane blocking: blocker non-fat dry milk (1g) in 1X Tris buffered saline (10ml, 1X TBS) + 0.1% Tween 20. After blocking (1 h), membrane washed with 1X Tris for 10 min to prepare for antibody.